



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Vortex Genie 2 T, 2700 rpm

**€559.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Vortex-Genie®, the universal small shaker with extra power!

High-speed vortexer.

Fast vortexing of problematic samples, easy to use. Circular motion, 6 mm orbit. Each including rubber plate (inside 2.5 cm Ø) and with 8.6 cm platform. Attachments easily exchangeable.

Order No. 7-0093 Vortex<sup>Genie®</sup>2-T

Press-to-Mix 1-60 sec., operation 1-60 min. with automatic stop after expiration.

Weight: 4 kg

Power supply: 240 V, 50 Hz

## Additional Information

No.	7-0093
Manufacturer (Brand)	neoLab
EAN	4058072027025
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	2700 rpm
Color	blue
sterile	No
Wide	122 mm
Height	165 mm
Length	165 mm
Depth	165 mm
Operation via	Knobs/Buttons
Weight	4 kg
Voltage	230 V
for ambient temperature	0 to 40 °C
Shaker properties	with infinitely variable speed control with timer with press-to-start function with continuous function including attachment (Further) attachments are available as accessories
Figure	2
Amplitude	1 mm
Number of speed levels	stepless
Attachment is suitable for	Reaction vessels. Test tubes. Centrifuge tubes (standard attachment)
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes
Movement	vibrating circulating
Speed MIN	600 rpm
Spare part	7-0510 7-0513 7-0514

Power	80 W
Maximum relative humidity at 20 °C	95.0 %
Model Shaker/ Rotator	neoLab Vortex Mixer Genie 2 T (7-0093)
Shaker type	Vibratory shaker/ vortexer

