




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Vortex Genie 2 Digital, 2850 rpm

€579.00
plus VAT &
Shipping

Product Images



Description

Vortex-Genie®, the universal small shaker with extra power!

High-speed vortexer.

Fast vortexing of problematic samples, easy to use. Circular motion, 6 mm orbit. Each including rubber plate (inside 2.5 cm Ø) and with 8.6 cm platform. Attachments easily exchangeable.

Order No. 7-0094 Vortex^{Genie®}2 -Digital

A digital display shows accuracy, reproducibility.

Press-to-Mix 1-90 sec., Operation 1-99 min.

Weight: 4 kg

Power supply: 240 V, 50 Hz

Additional Information

No.	7-0094
Manufacturer (Brand)	neoLab
EAN	4058072027032
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	2850 rpm
Color	blue
sterile	No
Wide	122 mm
Height	165 mm
Length	165 mm
Depth	165 mm
Operation via	Membrane keypad
Weight	4 kg
Voltage	230 V
for ambient temperature	0 to 40 °C
Shaker properties	programmable with infinitely variable speed control with timer with press-to-start function with speed monitoring with continuous function including attachment (Further) attachments are available as accessories
Amplitude	1 mm
Figure	3
Movement	vibrating circulating
Speed MIN	500 rpm
Number of speed levels	stepless
Power	80 W
Attachment is suitable for	Reaction vessels. Test tubes. Centrifuge tubes (standard attachment)
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes

Spare part	7-0510 7-0513 7-0514
Maximum relative humidity at 20 °C	95.0 %
Model Shaker/ Rotator	neoLab Vortex Mixer Genie 2 Digital (7-0094)
Shaker type	Vibratory shaker/ vortexer

