



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph Polymax 2040 (5° tilt angle) -- EU cable

**€2,922.00**  
**plus VAT & Shipping**

### Product Images



## Description

---

The loadable With a usable area of 390 x 340 mm, the Polymax 2040 is the optimal shaker for increased sample throughput. Since it can be equipped with a second shelf, the surface area can even be doubled.

## Additional Information

---

No.	BH-1749
Manufacturer (Brand)	Heidolph
VGKL number	524140005
EAN	4058072778989
Gross weight	15.62
Transport temperature	Room temperature
Net weight	13.06000 kg
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.52000 m
Packing height	0.33000 m
Packing depth	0.53000 m
Packaging volume	0.09100000 ccm
Customs tariff number	84798200

