



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Reax 2 -- EU cable

€1,235.00
plus VAT &
Shipping

Product Images



Description

The compact one for small volumes The Reax 2 with the optional adapter for test tubes is designed for simultaneous mixing of up to 20 vessels with diameters of 10–18 mm. This corresponds to the volume of classic test tubes or 15 ml Falcon tubes. Alternatively, you can clamp your complete rack into the universal adapter (width: 230 mm, height min./max.: 45 mm/160 mm) of the overhead shaker. Equipped with a solid metal frame, its construction stands for extreme resilience and reliability. It guarantees safety even during unattended continuous operation.

Additional Information

No.	BH-1737
Manufacturer (Brand)	Heidolph
EAN	4058072778873
Gross weight	5.4
Transport temperature	Room temperature
Net weight	4.80000 kg
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.23000 m
Packing height	0.26000 m
Packing depth	0.59000 m
Packaging volume	0.03500000 ccm
Customs tariff number	84798200

