



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Heidolph Reax 20 / 12 -- 1-16rpm -- EU cable

**€4,366.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

The loadable with a holding capacity of up to 12 bottles (image similar) The Reax 20 is an optimal shaker for processing and storing samples that tend to sediment. Applications range from environmental technology and paint processing to quality assurance in the materials industry. The maximum permissible intake weight is 30 kg, so that even samples with a high solids content are mixed properly. Since large forces are set in motion with this device, it is equipped with an emergency stop switch for your safety. This overhead shaker can be ordered with either a holding capacity for 4 , 8 or 12 bottles. Also suitable for mixing cylinders or wide neck bottles.

## Additional Information

---

No.	BH-1734
Manufacturer (Brand)	Heidolph
VGKL number	524120012
EAN	4058072778842
Gross weight	46
Transport temperature	Room temperature
Net weight	32.20000 kg
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.80000 m
Packing height	0.75000 m
Packing depth	1.20000 m
Packaging volume	0.72000000 cc
Customs tariff number	84798200

