



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## sunlab® 3D lab shaker with rubber mat 26 x 26 cm

**€430.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The Sunlab® 3D laboratory shaker is the optimal solution for various mixing tasks in different fields such as cell culture, microbiology or molecular biology. A fixed speed of 20 rpm as well as an inclination angle of 20° ensure a uniform tumbling motion and thus a homogeneous and efficient mixing of the samples. On the 26 x 26 cm platform of the Sunlab® 3D laboratory shaker, a rubber mat with elastic nubs ensures a secure and firm hold of sample vessels such as reaction tubes, vials, dishes or plates.

Application:

- Staining and decolorizing gels,
- mixing protein solutions,
- hybridizing DNA, etc.

## Additional Information

No.	D-8040
Manufacturer (Brand)	sunlab
EAN	4058072091729
Transport temperature	Room temperature
Speed MAX	20 rpm
Color	orange
Wide	260 mm
Height	190 mm
Length	260 mm
Depth	260 mm
Shaker properties	with permanent function including rubber mat(s)
Operation via	Knobs/Buttons
Movement	circling staggering
Width tray	260 mm
Speed MIN	20 rpm
Weight	1.8 kg
Power	14 W
Length shelf	260 mm
Model Shaker/ Rotator	Sunlab 3D Orbital Shaker (D-8040)
Inclination angle	20°
Protection class	IP 21
Voltage	230 V
Sunlab® Model number	Shaker SUI040
Shaker type	Orbital/Circular Shaker Rocker Shaker
loadable up to	0.8 kg
for ambient temperature	5 to 45 °C

