




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Rocker shaker, adjustable inclination

€750.40
plus VAT &
Shipping

Product Images



Description

Rocker shaker with variable rocking frequency

The large platform rocks along a central axis. With non-slip rubber mat for Petri dishes, microtiter plates or small vessels.

Ideally suited for Western blots.

The tilt angle is 4-fold adjustable (3, 6, 9, 12°)

Technical data

Weight: 3.8 kg

Connection: 220 Volt 50 Hz

Additional Information

No.	E-1113
Manufacturer (Brand)	neoLab
EAN	4058072062941
Transport temperature	Room temperature
Speed MAX	100 rpm
Color	White
Wide	350 mm
Height	110 mm
Length	520 mm
Depth	520 mm
Shaker properties	with infinitely variable speed control with speed monitoring with continuous function including rubber mat(s)
Number of speed levels	stepless
Operation via	Knobs/Buttons
Movement	teetering
Width tray	340 mm
Speed MIN	10 rpm
Weight	3.8 kg
Power	7 W
Length shelf	410 mm
Model Shaker/ Rotator	neoLab Rocker Shaker (E-1113)
MAX tilt angle	12 °
MIN tilt angle	3 °
Voltage	230 V
Shaker type	Rocker shaker
loadable up to	3.8 kg
for ambient temperature	2 to 50 °C

