



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN analytical balance ABJ, 220 g

€1,500.00
plus VAT &
Shipping



Description

- ABJ-NM only: Internal automatic calibration/adjustment when temperature changes $> 2\text{ }^{\circ}\text{C}$ and time-controlled every 4 hours, guarantees high accuracy and makes calibration/adjustment independent of location
- ABS-N only: CAL adjustment program for setting accuracy using an external test weight
- Dosing aid: High stability mode and other filter settings selectable
- Convenient formulation/documentation with combined tare/print function. In addition, the mixture ingredients of the formulation are automatically numbered and printed with number/weight value
- Scale identification number: 4-digit, freely programmable, is printed in the calibration record
- Automatic data output to PC/printer after each balance stability
- ABJ-NM has OIML certificate

Single-cell advanced technology:

- Fully automatic production of the load cell from one piece
- Stable temperature behavior
- Short settling time: Stable weighing values within approx. 3 sec under laboratory conditions
- High mechanical robustness
- High corner load safety
- Dust cover included in delivery

Additional Information

No.	KP-0758
Manufacturer (Brand)	KERN
VGKL number	105835220
EAN	4058072126568
Transport temperature	Room temperature
Color	light greylgrey
Wide	214 mm
Height	340 mm
Depth	340 mm
Readability	0,1 mg
Operation via	Membrane keypad
Weight	9.22 kg
Reproducibility	0,2 mg
Scales features	with switchable weighing units
Shape weighing plate	around
Material weighing plate	Stainless steel
Weighing range MAX	220 g
Manufacturer	CORE
Type scale	Analytical balance

