



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN analytical balance ALJ, 510 g

€2,250.00
plus VAT &
Shipping

Product Images



Description

- KERN ALJ 200-5DA: High precision semi-micro analytical balance
- KERN ALJ-A03: Ionizer for neutralizing electrostatic charge for fixed installation in the analytical balance. This makes handling particularly convenient, as a separate device is no longer necessary. Simply press a button to switch on the ionization blower. Suitable for all models of this series
- aLJ-A/-AM only: Internal automatic adjustment for temperature changes ± 1.2 °C and time-controlled according to defined intervals, guarantees high accuracy and makes the unit location-independent
- aLS-A only: CAL adjustment program for setting accuracy, using an external test weight
- Ergonomically optimized control panel for left- and right-handed operation
- Large glass draft shield with 3 sliding doors for convenient access to the weighing product
- Compact dimensions, advantageous where space is limited
- Dust cover included in delivery

Additional Information

No.	KP-0575
Manufacturer (Brand)	KERN
VGKL number	105838504
EAN	4058072126094
Transport temperature	Room temperature
Color	light grey grey
Wide	210 mm
Height	340 mm
Depth	340 mm
Readability	0,1 mg
Operation via	Membrane keypad
Weight	10.7 kg
Reproducibility	0,2 mg
Type interface	RS232
Scales features	with switchable weighing units with internal calibration with draft shield with interface with GLP functions
Energy supply	Battery Mains
Shape weighing plate	around
Material weighing plate	Stainless steel
Calibration type	at the touch of a button
Weighing range MAX	510 g
Manufacturer	CORE
Type scale	Analytical balance
Weighing units	ct lb oz

