



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN analytical balance ALJ, 160 g

**€1,460.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- KERN ALJ 200-5DA: High precision semi-micro analytical balance
- KERN ALJ-A03: Ionizer for neutralizing electrostatic charge for fixed installation in the analytical balance. This makes handling particularly convenient, as a separate device is no longer necessary. Simply press a button to switch on the ionization blower. Suitable for all models of this series
- aLJ-A/-AM only: Internal automatic adjustment for temperature changes  $\pm 1.2$  °C and time-controlled according to defined intervals, guarantees high accuracy and makes the unit location-independent
- aLS-A only: CAL adjustment program for setting accuracy, using an external test weight
- Ergonomically optimized control panel for left- and right-handed operation
- Large glass draft shield with 3 sliding doors for convenient access to the weighing product
- Compact dimensions, advantageous where space is limited
- Dust cover included in delivery

## Additional Information

---

No.	KP-0619
Manufacturer (Brand)	KERN
VGKL number	105838165
EAN	4058072126247
Transport temperature	Room temperature
Color	light grey grey
Wide	215 mm
Height	345 mm
Depth	345 mm
Readability	0,1 mg
Operation via	Membrane keypad
Weight	9.2 kg
Reproducibility	0,1 mg
Type interface	RS232
Scales features	with switchable weighing units with internal calibration with draft shield with interface with GLP functions
Energy supply	Battery Mains
Shape weighing plate	around
Manufacturer	CORE
Material weighing plate	Stainless steel
Calibration type	at the touch of a button
Type scale	Analytical balance
Weighing range MAX	160 g
Weighing units	ct lb oz

