



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Sartorius Secura® semi-microbalance, 120 | 220 g, weighing plate 80 mm Ø, fully automatic calibration and adjustment

**€8,400.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

**Secura® scales** with integrated and intelligent safety systems.

Secura® scales provide you with highly accurate weighing results and ease of use, as well as integrated safety systems for all-round reliability – simply weigh safely.

Real-time operator guidance when leveling the scale, alarm indications and interactive assistance support leveling, thus you avoid inaccurate results due to an improperly leveled scale.

The internal fully automatic calibration/adjustment isoCAL for highest possible accuracy can be performed at any time. Cal Audit Trail enables easy and traceable documentation of all calibration and adjustment procedures. All weighing data are traceable thanks to secure and clear documentation with sample and batch ID numbers, GLP/GMP-compliant printouts. Standard password protection prevents tampering and unauthorized access. Easy record keeping with the integrated PC Direct function and Plug & Work technology for connection to a Sartorius laboratory printer.

## Additional Information

No.	ST-8801
Manufacturer (Brand)	Sartorius
old neoLab article no.	5-2040
EAN	4058072371371
Transport temperature	Room temperature
Color	gray white
Material	Acrylonitrile butadiene styrene (ABS)
Wide	218 mm
Height	316 mm
Length	376 mm
Depth	376 mm
Scales features	with switchable weighing units with internal calibration with draft shield with interface including batteries/accumulator
Readability	0.0100 g
Operation via	Touchscreen
DM Weighing plate	80 mm
Energy supply	Power grid
Shape weighing plate	around
Weight	7.8 kg
Material weighing plate	Stainless steel
Reproducibility	0.03 g
Calibration type	automatic
Type interface	USB mini
Type scale	Analytical balance
Weighing range MAX	220 g
Weighing units	g ct mommes kg lb ozt oz baht grain
for ambient temperature	10 to 30°C

