



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN spring scale 281, 300 g

€105.00
plus VAT &
Shipping

Product Images



Description

- **Scale tube made of aluminum:** robust, durable, rustproof
- **Gram/kilogram scale:** Measurement result displayed in grams or kilograms instead of N
- **Compressive force measurement:** possible by means of an optional pressure set
- **Retaining bracket:** thanks to the rotatable retaining bracket, the scale can always be optimally aligned in the direction of view
- **High precision:** backlash-free spring bearing with integrated tare screw for high-precision adjustment
- Fatigue-free stainless steel spring
- Clip eyelet with free rotation of the lower suspension by 360° for models with ≤ 1 kg
- **High-quality workmanship:** abrasion-resistant, color-anodized precision scale with high resolution for the finest reading of the measurement result

Additional Information

No.	KP-0918
Manufacturer (Brand)	KERN
EAN	4058072126643
Transport temperature	Room temperature
Color	light greylgrey
Wide	220 mm
Height	12 mm
Depth	12 mm
Operation via	read only
Weight	47 kg
Scales features	with switchable weighing units
Weighing range MAX	300 g
Manufacturer	CORE
Type scale	Spring scale

