



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN spring scale 281, 60 g

€92.00
plus VAT &
Shipping

Product Images



Description

- **Scale tube made of aluminum:** robust, durable, rustproof
- **Gram/kilogram scale:** Measurement result displayed in grams or kilograms instead of N
- **Compressive force measurement:** possible by means of an optional pressure set
- **Retaining bracket:** thanks to the rotatable retaining bracket, the scale can always be optimally aligned in the direction of view
- **High precision:** backlash-free spring bearing with integrated tare screw for high-precision adjustment
- Fatigue-free stainless steel spring
- Clip eyelet with free rotation of the lower suspension by 360° for models with ≤ 1 kg
- **High-quality workmanship:** abrasion-resistant, color-anodized precision scale with high resolution for the finest reading of the measurement result

Additional Information

| | |
|-----------------------|------------------|
| No. | KP-3867 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072135669 |
| Transport temperature | Room temperature |
| Color | light greylgrey |
| Wide | 220 mm |
| Height | 12 mm |
| Depth | 12 mm |
| Operation via | Membrane keypad |
| Weight | 39 kg |
| Manufacturer | CORE |
| Type scale | Spring scale |
| Weighing range MAX | 60 g |

