



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN spring scale 285, 50 kg

€255.00
plus VAT &
Shipping

Product Images



Description

- **Scale tube made of aluminum:** robust, durable, rustproof
- **Gram/kilogram scale:** Measurement result displayed in grams or kilograms instead of N
- **Compressive force measurement:** possible by means of an optional pressure set
- **Drag pointer and handle:** standard on all Series 285 models
- **Retaining bracket:** thanks to the rotatable retaining bracket, the scale can always be optimally aligned in the direction of view
- **High precision:** backlash-free spring bearing with integrated tare screw for high-precision adjustment
- Fatigue-free stainless steel spring
- Clip eyelet with free rotation of the lower suspension by 360° for models with ≤ 1 kg
- **High-quality workmanship:** abrasion-resistant, color-anodized precision scale with high resolution for the finest reading of the measurement result

Additional Information

No.	KP-1559
Manufacturer (Brand)	KERN
EAN	4058072128227
Transport temperature	Room temperature
Color	light greylgrey
Wide	370 mm
Height	32 mm
Depth	32 mm
Operation via	read only
Weight	570 kg
Scales features	with switchable weighing units
Weighing range MAX	50 kg
Manufacturer	CORE
Type scale	Spring scale

