



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN hanging scale HCB, 500 kg

**€215.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With the TÜV mark of conformity, the scales meet the requirements of EN 13155 (loose load receptors) and EN 61010-1 (electrical safety) Ideal for fast checkweighing in incoming and outgoing goods Indispensable also in the private sector for determining the weight of fish, game, fruit, bicycle parts, suitcases, etc. Hold function: For convenient reading of the weight value, the display can be "frozen" in various ways. Either automatically when the scale stops weighing or manually by pressing the hold button Peak load display (peak hold), measuring frequency 5 Hz HCB with  $\geq 0.5$  t: Reinforced version for high loads. Carabiner and hook not standard

## Additional Information

---

No.	KP-4334
Manufacturer (Brand)	KERN
EAN	4058072139629
Transport temperature	Room temperature
Color	light gray dark gray
Readability	1 kg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Type scale	Suspension scale
Weighing range MAX	500 kg

