



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN hanging scale HCB, 100 kg

**€110.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- With the TÜV mark of conformity, the scales meet the requirements of EN 13155 (loose load receptors) and EN 61010-1 (electrical safety) standards
- Ideal for fast checkweighing in incoming and outgoing goods departments
- Indispensable also in the private sector for determining the weight of fish, game, fruit, bicycle parts, suitcases, etc.
- Hold function: For convenient reading of the weight value, the display can be "frozen" in various ways. Either automatically when the scale stops weighing or manually by pressing the Hold button
- Peak load display (peak hold), measuring frequency 5 Hz
- HCB with  $\geq 0.5$  t: Reinforced version for high loads. Carabiner and hook not standard

## Additional Information

---

No.	KP-4335
Manufacturer (Brand)	KERN
VGKL number	645875100
EAN	4058072139636
Transport temperature	Room temperature
Color	light gray dark gray
Readability	200 g
Operation via	Membrane keypad
Energy supply	Power grid
Weighing range MAX	100 kg
Manufacturer	CORE
Type scale	Suspension scale

