



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN hanging scale CH , 15 kg

€75.00
plus VAT &
Shipping

Product Images



Description

- With the TÜV mark of conformity, the scales meet the requirements of EN 13155 (loose load receptors) and EN 61010-1 (electrical safety) standards
- Ideal for fast checkweighing in incoming and outgoing goods departments
- Indispensable also in the private sector for weight determination of fish, game, fruit, bicycle parts, suitcases, etc.
- Hold function: For convenient reading of the weight value, the display can be "frozen" in various ways. Either automatically when the scale stops weighing or manually by pressing the hold button
- Measuring tape, pull-out, measuring length approx. 100 cm
- Peak load display (peak hold)
- Ready to start: Batteries included, 9 V block, operating time up to 20 h

Additional Information

No.	KP-1998
Manufacturer (Brand)	KERN
VGKL number	645871015
EAN	4058072129361
Transport temperature	Room temperature
Color	light gray dark gray
Readability	20 g
Operation via	Membrane keypad
Energy supply	Power grid
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Suspension scale

