



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN hanging scale HDB, 10 kg

**€65.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- With the TÜV mark of conformity, the scales meet the requirements of EN 13155 (loose load receptors) and EN 61010-1 (electrical safety)
- NEW: KERN HDB now also available as high capacity XL version. The larger design creates space for a large handle, which allows even better handling especially with heavier loads
- Ideal for fast checkweighing in incoming and outgoing goods departments
- Indispensable also in the private sector for determining the weight of fish, game, fruit, bicycle parts, suitcases, etc.
- Hold function: For convenient reading of the weight value, the display can be "frozen" in various ways. Either automatically when the scale stops weighing or manually by pressing the Hold key
- Peak load display (peak hold)
- Hook, steel, foldable

## Additional Information

---

No.	KP-4355
Manufacturer (Brand)	KERN
EAN	4058072139834
Transport temperature	Room temperature
Color	light gray dark gray
Readability	10 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Type scale	Suspension scale
Weighing range MAX	10 kg

