



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Luggage scale HGA, 50 kg

**Price on
request
plus VAT &
Shipping**



Description

- This ergonomically optimized luggage scale is ideal for quick and precise checking of suitcases, bags, sports luggage, etc.. The high maximum load ensures sufficient safety reserve and is therefore also suitable for weighing other hanging loads
- With the TÜV test mark, the scale complies with the requirements of the standard: EN 13155 (loose load receptors/breakage safety) and EN 61010-1 (electrical safety)
- Ergonomic design for comfortable handling
- Sturdy nylon webbing strap for secure fastening of the weighing product
- Simple and convenient 1-key operation
- Convenient unit changeover at the touch of a button
- High mobility: battery operation and compact, lightweight design make it suitable for use at multiple locations
- Automatic data hold function: when the scale stops weighing, the weight display is "frozen" for 12 sec. This leaves enough time to read the weight value at rest
- Ready to start: Batteries included, 3V CR2032, operating time up to 20 h, AUTO-OFF function for battery saving. Battery level indicator in the display
- Note: The models are supplied in sets of 5 only. I.e. the prices given in the table refer to a delivery of 5 pieces each. Individual delivery not possible. The calibration prices given refer to the calibration of a single scale

Additional Information

No.	KP-4397
Manufacturer (Brand)	KERN
EAN	4058072140205
Transport temperature	Room temperature
Color	light gray dark gray
Readability	10 g
Operation via	Membrane keypad
Energy supply	Power grid
Weighing range MAX	50 kg
Manufacturer	CORE
Type scale	Suitcase scale

