



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN pallet scale TGD 1500 kg

**Price on
request
plus VAT &
Shipping**

Product Images



Description

High mobility: thanks to battery operation (optional) and compact, flat design suitable for use at several locations (production, warehouse, shipping, etc.) Load capacity: steel, powder-coated, 4 silicone-coated steel load cells, dust and splash protection IP67. Also available as a component without evaluation unit, evaluation unit Two castors and handle for convenient transport of the scale Weighing with tolerance range (checkweighing): A visual and acoustic signal supports fast checking of palletized goods Totalizing weight values and counting parts Protective work hood included Did you know? Our floor scales are delivered in a robust wooden box. This protects the high quality weighing technology from environmental influences and stress during transport. KERN - always one idea ahead Shipping by forwarding agency. Please ask for dimensions, gross weight, shipping costs

Additional Information

No.	KP-5568
Manufacturer (Brand)	KERN
EAN	4058072150730
Transport temperature	Room temperature
Color	light gray dark gray
Readability	500 g
Operation via	Membrane keypad
Energy supply	Power grid
Weighing range MAX	1 t
Manufacturer	CORE
Type scale	Pallet scale

