



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**KERN Mechanical weight and force
measurement with quality spring for long life
289-100,**

**€81.00
plus VAT &
Shipping**

Product Images



Description

- Optimized price/performance ratio thanks to transparent plastic housing, ideal for schools and educational institutions
- Newton scale: The KERN 281 series features a result display in newtons instead of grams, especially for measurements of tensile forces
- High quality:
 - - Backlash-free spring bearing with high-precision adjustment
 - - Fatigue-free Inox spring
- Abrasion-free, colored precision scale with high resolution
- Thanks to the rotatable inner tube, the scale can always be read optimally
- The clamp supplied as standard can be conveniently exchanged for another suspension device, thus individually adaptable to the weighing product
- A large number of other spring scales with gram or Newton graduation for tension and compression measurements and specific accessories can be found at www.sauter.eu

Additional Information

No.	KP-5657
Manufacturer (Brand)	KERN
EAN	4058072151317
Transport temperature	Room temperature
Color	light gray
Readability	0,01 N
Weighing range MAX	1 N
Manufacturer	CORE
Type scale	Force gauge
Weighing units	Newton

