



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**KERN force test stand - electronic, one column
Max 500 N: Vertical <-> 300 mm**

**€2,480.00
plus VAT &
Shipping**

Product Images



Description

Motorized test stand for tensile and compressive force tests Table design for easy operation Robust design for durable use Easy-to-operate stop switch Individually adjustable end points of travel Automatic or manual travel mode Suitable for all SAUTER force gauges up to 500 N (e.g. , not included)

Additional Information

No.	KP-5553
Manufacturer (Brand)	KERN
EAN	4058072150587
Transport temperature	Room temperature
Color	light gray
Weighing range MAX	500 N
Manufacturer	CORE
Type scale	Force gauge
Weighing units	Newton

