



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN elec. Force gauge (internal) Max 1 kN: d=0.5 N

€275.00
plus VAT &
Shipping

Product Images



Description

Reversible display: automatic direction recognition Safe operability due to ergonomic housing design Peak hold function for recording the peak value or track function for continuous measurement display Selectable units: N, lb, kg, oz Auto-Power-Off Standard attachments: as shown, extension rod: 90 mm Can be attached to all SAUTER test stands With one of the two optional attachments for tensile stress testing, the SAUTER FK becomes a tensiometer for testing material stress properties of cables, wires, threads, cords, etc. (see accessories)

Additional Information

No.	KP-4243
Manufacturer (Brand)	KERN
EAN	4058072138769
Transport temperature	Room temperature
Color	light gray
Readability	0,5 N
Weighing range MAX	1 N
Manufacturer	CORE
Type scale	Force gauge
Weighing units	Newton

