



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN elec. Force gauge (internal) Max 500 N: d=0.2 N

**€275.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Reversible display: automatic direction recognition Safe operability due to ergonomic housing design Peak hold function for recording the peak value or track function for continuous measurement display Selectable units: N, lb, kg, oz Auto-Power-Off Standard attachments: as shown, extension rod: 90 mm Can be attached to all SAUTER test stands With one of the two optional attachments for tensile stress testing, the SAUTER FK becomes a tensiometer for testing material stress properties of cables, wires, threads, cords, etc. (see accessories)

## Additional Information

---

No.	KP-4247
Manufacturer (Brand)	KERN
EAN	4058072138806
Transport temperature	Room temperature
Color	light gray
Readability	0,2 N
Weighing range MAX	500 N
Manufacturer	CORE
Type scale	Force gauge
Weighing units	Newton

