



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN mech. force gauge Max 20 N: d=0.1 N

**€230.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Dual scale: Newton and kg can be read in parallel Rotatable reading unit for easy adjustment of the instrument Peak-hold function by drag pointer Mountable on all manual test stands Zero setting of the display at the push of a button Supplied in a robust carrying case Standard attachments: as shown, extension rod: 90 mm

## Additional Information

---

No.	KP-4165
Manufacturer (Brand)	KERN
EAN	4058072138158
Transport temperature	Room temperature
Color	light gray
Readability	0,1 N
Weighing range MAX	20 N
Manufacturer	CORE
Type scale	Force gauge
Weighing units	Newton

