




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN mech. force gauge Max 500 N: d=2.5 N

€230.00
plus VAT &
Shipping

Product Images



Description

Dual scale: Newton and kg can be read in parallel Rotatable reading unit for easy adjustment of the instrument Peak-hold function by drag pointer Mountable on all manual test stands Zero setting of the display at the push of a button Supplied in a robust carrying case Standard attachments: as shown, extension rod: 90 mm

Additional Information

No.	KP-4170
Manufacturer (Brand)	KERN
EAN	4058072138202
Transport temperature	Room temperature
Color	light gray
Readability	2,5 N
Manufacturer	CORE
Type scale	Force gauge
Weighing range MAX	500 N
Weighing units	Newton

