



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN Baby Scale MBC, 20 kg

€410.00
plus VAT &
Shipping

Product Images



Description

- Calibration class III (calibration optional)
- Medical approval according to 93/42/EEC
- Popular in pediatrics thanks to the large bowl, which is ideal for infants with its optimal ergonomic curvature
- Special function for neonatology wards: Differential function for calculating weight gain before and after feeding
- Robust ABS plastic. Hygienic and easy to clean
- Safe and slip-free positioning due to rubber feet
- Hold function: In the case of restless patients, a stable weight value is determined and frozen by averaging. This leaves enough time to devote to the small patient and then read off the weight value in peace
- Battery and mains operation possible, battery operation optional
- In the EU, calibration is required by law for scales intended for use as medical devices, please order at the same time. For this we need installation site and postal code.
- NEW: MBC-E, identical to MBC, but without power adapter, especially suitable for mobile use, e.g. for midwives
- KERN MBC: battery and mains operation possible, battery operation optional
- MBC-E: battery operation possible, mains and battery operation optional

Additional Information

No.	KP-4692
Manufacturer (Brand)	KERN
EAN	4058072142926
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	10 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Plastic
Weighing range MAX	20 kg

