



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN personal scale MPB, 300 kg

€230.00
plus VAT &
Shipping

Product Images



Description

- Sturdy steel weighing plate with non-slip and abrasion-resistant surface
- Hygienic and easy to clean
- Safe and slip-free positioning due to rubber feet
- Dynamic weighing function to determine a stable weighing value, ideal for restless patients
- Hold function: For restless patients, a stable weight value is determined and "frozen" by averaging. This leaves enough time to devote to the patient and then read off the weight value in peace
- Simple and convenient 4-key operation
- BMI function for determination of underweight/normal weight/overweight
- Battery and mains operation possible
- Wall bracket for wall mounting of the indicator as standard, only for model without stand

Additional Information

No.	KP-4717
Manufacturer (Brand)	KERN
EAN	4058072153977
Transport temperature	Room temperature
Color	light gray dark gray
Readability	100 g
Operation via	Membrane keypad
Energy supply	Power grid
Weighing range MAX	300 kg
Manufacturer	CORE
Type scale	Personal scale

