



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN personal scale MPC, 250 kg

**€400.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Version with = 300 kg and larger weighing plate! Particularly suitable for weighing overweight patients in the obesity range
- Calibration class III (calibration optional)
- Medical approval according to 93/42/EWG
- The bestseller for medical practices, pharmacies and hospitals
- Robust design for daily use in professional environments
- Ergonomically optimized design
- Hygienic and easy to clean
- Safe and slip-free set-up due to height-adjustable rubber feet
- Bubble level for exact leveling of the scale as standard equipment
- Sturdy handle for transporting the scale
- Evaluation unit can be positioned in many ways: e.g. free-standing or screwed to the wall. Wall bracket for mounting the indicator on the wall as standard equipment
- Large tread surface with non-slip and abrasion-resistant plastic surface, for optimum patient stability
- Hold function: If the patient is restless, a stable weight value is determined by averaging and "frozen". This leaves enough time to attend to the patient and then read off the weight value in peace
- BMI function for determining underweight/normal weight/overweight
- Mother-child function: the weight of the mother determined by weighing is set to zero at the touch of a button. The mother and child are then weighed together. The child's weight is now shown directly on the display. This is particularly convenient for infants, who can be weighed safely and securely on their mother's arm
- Fine display: by pressing a button, the readability can be increased by 1 digit for 5 seconds
- Battery and mains operation possible, battery operation optional
- In the EU, calibration is required by law for scales intended for use as medical devices, please order at the same time. For this we need installation site and zip code.

## Additional Information

---

No.	KP-4721
Manufacturer (Brand)	KERN
EAN	4058072143152
Transport temperature	Room temperature
Color	light gray dark gray
Readability	100 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Plastic
Type scale	Personal scale
Weighing range MAX	250 kg

