



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN personal scale MPE, 250 kg

€615.00
plus VAT &
Shipping

Product Images



Description

- Calibration class III
- Medical approval according to 93/42/EWG
- Personal scale KERN MPE-E identical in construction with MPE, but without mains adapter, only battery operation
- Practical second display on the back of the scale, therefore easy to read by the patient and the doctor or nursing staff
- The MPE-HM model has an extendable size measuring rod elegantly integrated in the stand, with a measuring range of 88 - 205 cm
- The ergonomically optimized design and the robust construction predestine this series for daily use in professional environments
- Safe and slip-free set-up due to height-adjustable rubber feet
- Bubble level for exact leveling of the scale as standard equipment
- Large tread surface with non-slip and abrasion-resistant plastic surface for optimum patient stability
- Hold function
- Mother-child function: the weight of the mother determined by weighing is set to zero at the touch of a button. The mother and child are then weighed together. The child's weight is now shown directly on the display. This is particularly convenient for infants, who can be weighed safely and securely on their mother's arm
- BMI function for determining normal weight/overweight
- Fine display: by pressing a button, the readability can be increased by 1 digit for 5 seconds
- Convenient to transport thanks to the two integrated transport wheels
- Battery and mains operation possible, battery operation optional

Additional Information

No.	KP-4727
Manufacturer (Brand)	KERN
EAN	4058072143206
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	100 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Plastic
Weighing range MAX	250 kg
Manufacturer	CORE
Type scale	Personal scale

