



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN personal scale MPE, 250 kg

€665.00
plus VAT &
Shipping

Product Images





Description

- Calibration class III (calibration optional)
- Medical approval according to 93/42/EEC
- Practical second display on the back of the scale, thus conveniently readable by the patient and the doctor or nursing staff
- The MPE-HM model has an extendable size measuring rod elegantly integrated in the stand, with a measuring range of 88 - 205 cm
- The ergonomically optimized design and the robust construction predestine this series for daily use in professional environments
- Hygienic and easy to clean
- Safe and slip-free set-up due to height-adjustable rubber feet
- Level for exact levelling of the scale as standard
- Large tread surface with non-slip plastic surface
- Vibration-free weighing: Vibrations are filtered out to determine a stable weighing value, ideal for restless standing patients
- Hold function: For restless patients, a stable weight value is determined and "frozen" by averaging. This leaves enough time to devote to the patient and then read off the weight value in peace
- Mother-child function: the weight of the mother determined by weighing is set to zero at the touch of a button. The mother and child are then weighed together. The child's weight is now shown directly on the display. Particularly convenient for infants, who can be weighed safely and securely on the mother's arm
- BMI function for determining underweight/normal weight/overweight
- Convenient to transport thanks to the two integrated transport wheels
- Fine display: by pressing the button, the readability can be increased by 1 digit for 5 seconds
- Battery and mains operation possible, battery operation optional
- In the EU, calibration is required by law for scales intended for use as medical devices, please order at the same time. For this we need installation site and zip code.

Additional Information

No.	KP-4729
Manufacturer (Brand)	KERN
EAN	4058072143220
Transport temperature	Room temperature
Color	light gray dark gray
Readability	100 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Plastic
Weighing range MAX	250 kg
Manufacturer	CORE
Type scale	Personal scale

