



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN personal scale MPS, 200 kg

**€730.00**  
**plus VAT &**  
**Shipping**



## Description

---

- Calibration class III (calibration optional)
  - Medical approval according to 93/42/EEC
  - Thanks to its compact dimensions and good price/performance ratio, this series is ideal for treatment rooms in doctors' surgeries and hospitals
  - Ergonomic evaluation device with large keypad and high-contrast LCD display for convenient input and reading of patient values, such as height for calculating BMI
  - Sturdy steel weighing plate with non-slip and abrasion-resistant surface
  - Hygienic and easy to clean
  - Safe and slip-free set-up due to height-adjustable rubber feet
  - Level indicator for exact leveling of the scale as standard equipment
  - Dynamic weighing function for determining a stable weight value, ideal for patients who are standing unsteadily
  - Hold function: For restless patients, a stable weight value is determined and "frozen" by averaging. This leaves enough time to devote to the patient and then read off the weight value in peace
  - Intuitive operation via keypad
  - PRE-TARE function for 5 weight values
  - BMI function for determination of underweight/normal weight/overweight
  - Battery and mains operation possible
  - Wall bracket for wall mounting of the indicator as standard equipment
  - Battery and mains operation possible, battery operation optional
- In the EU, calibration is required by law for scales intended for use as a medical device, please order at the same time. For this we need the place of installation and zip code.

## Additional Information

---

No.	KP-4736
Manufacturer (Brand)	KERN
EAN	4058072143282
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	100 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Type scale	Personal scale
Weighing range MAX	200 kg

