



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN chair scale MCB, 300 kg

€1,400.00
plus VAT &
Shipping

Product Images



Description

- Calibration class III (calibration optional)
- Medical approval according to 93/42/EEC
- Mobile version with parking brakes on all wheels
- Due to their large cross-section, the castors make it easy to negotiate door thresholds, edges and gaps in passenger elevators
- Hold function: In the case of restless patients, a stable weight value is determined and frozen by averaging. This leaves enough time to attend to the patient and then read off the weight value in peace and quiet
- BMI function for determining underweight/normal weight/overweight
- This chair scale is the ideal measuring instrument for old people's homes, rehabilitation centers and clinics with a high volume of patients, including weighing in the obese range up to 300 kg
- With its four wheels and low dead weight, this chair scale ensures maximum mobility when moving to the patient. This results in increased time efficiency for clinical staff and improved safety for patients, who can be weighed in their familiar surroundings
- For frail patients, the comfortable chair shape provides secure support during weighing
- This chair scale is designed according to industry requirements, ensuring high long-term stability
- Two fold-down armrests and footrests make it easy to take a seat. Ideal for overweight patients or for barrier-free transfer, e.g. from bed to chair scale
- Battery and mains operation possible,

In the EU, calibration is required by law for scales intended for use as medical devices, please order at the same time. For this we need installation site and zip code. battery operation optional

Additional Information

No.	KP-4702
Manufacturer (Brand)	KERN
EAN	4058072143008
Transport temperature	Room temperature
Color	light gray dark gray
Readability	100 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Plastic
Weighing range MAX	300 kg
Manufacturer	CORE
Type scale	Chair scale

