




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN platform scale DE , 60 kg | 150 kg

**€225.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- High mobility: thanks to battery operation (optional), compact, flat design and low weight, suitable for use at several locations (production, warehouse, shipping, etc.)
- Evaluation unit can be positioned in a variety of ways: e.g. free-standing or screwed to the wall (optional)
- Evaluation unit: plastic, dust and splash water protection IP65
- Weighing plate: stainless steel, base steel painted
- PRE-TARE function for manual preweighing of a known container weight, useful for fill quantity checks
- With the recipe function, different ingredients of a mixture can be weighed. The total weight of all ingredients can be called up for checking purposes
- Freely programmable weighing unit, e.g. display directly in thread length g/m, paper weight g/m<sup>2</sup>, etc.
- Protective hood included in scope of delivery

## Additional Information

---

No.	KP-3985
Manufacturer (Brand)	KERN
VGKL number	105836037
EAN	4058072136628
Transport temperature	Room temperature
Color	light gray dark gray
Readability	20 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Platform scale
Weighing range MAX	150 kg

