



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN platform scale DE, 6 kg, 15 kg

€380.00
plus VAT &
Shipping

Product Images



Description

High mobility: Thanks to battery operation (optional), compact, flat design and low weight, suitable for use at several locations (production, warehouse, shipping, etc.) Evaluation unit can be positioned in many ways: e.g., free-standing or screwed to the wall (optional) Evaluation unit: Plastic, dust and splash water protection IP65 Weighing plate: Stainless steel, base steel painted PRE-TARE function for manual pre-drawing of a known container weight, useful for fill quantity checks The recipe function allows various components of a mixture to be weighed. The total weight of all ingredients can be called up for checking Freely programmable weighing unit, e.g. display directly in yarn length g/m, paper weight g/m², etc. Protective hood included in scope of delivery

Additional Information

No.	KP-3990
Manufacturer (Brand)	KERN
VGKL number	105836047
EAN	4058072136673
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	200 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	15 kg
Manufacturer	CORE
Type scale	Platform scale

