



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN platform scale DE, 15 kg, 35 kg

**€375.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

High mobility: Thanks to battery operation (optional), compact, flat design and low weight, suitable for use at several locations (production, warehouse, shipping, etc.) Evaluation unit can be positioned in many ways: e.g., free-standing or screwed to the wall (optional) Evaluation unit: Plastic, dust and splash water protection IP65 Weighing plate: Stainless steel, base steel painted PRE-TARE function for manual pre-drawing of a known container weight, useful for fill quantity checks The recipe function allows various components of a mixture to be weighed. The total weight of all ingredients can be called up for checking Freely programmable weighing unit, e.g. display directly in yarn length g/m, paper weight g/m<sup>2</sup>, etc. Protective cover included in scope of delivery

## Additional Information

---

No.	KP-3996
Manufacturer (Brand)	KERN
VGKL number	105836048
EAN	4058072136734
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	500 mg
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	35 kg
Manufacturer	CORE
Type scale	Platform scale

