



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN platform scale DE , 3 kg | 6 kg

€210.00
plus VAT &
Shipping

Product Images



Description

- High mobility: thanks to battery operation (optional), compact, flat design and low weight, suitable for use at several locations (production, warehouse, shipping, etc.)
- Evaluation unit can be positioned in a variety of ways: e.g. free-standing or screwed to the wall (optional)
- Evaluation unit: plastic, dust and splash water protection IP65
- Weighing plate: stainless steel, base steel painted
- PRE-TARE function for manual preweighing of a known container weight, useful for fill quantity checks
- With the recipe function, different ingredients of a mixture can be weighed. The total weight of all ingredients can be called up for checking purposes
- Freely programmable weighing unit, e.g. display directly in thread length g/m, paper weight g/m², etc.
- Protective cover included in scope of delivery

Additional Information

No.	KP-4005
Manufacturer (Brand)	KERN
VGKL number	105836031
EAN	4058072136826
Transport temperature	Room temperature
Color	light gray dark gray
Readability	1 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	6 kg
Manufacturer	CORE
Type scale	Platform scale

