



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN platform scale DS , 150 kg

€1,140.00
plus VAT &
Shipping

Product Images



Description

- High-load precision balance, ideal for large-volume or heavy samples that must be weighed with high accuracy
- Step-by-step operation via yes/no dialog on display
- Numerical tare weight for known container weight. Useful for fill quantity checks
- Accurate counting: Automatic reference sample updating gradually improves the average value of the part weight
- Freely programmable weighing unit, e.g. display directly in thread length g/m, paper weight g/m², etc.
- Protective work cover included in scope of delivery

Additional Information

No.	KP-4013
Manufacturer (Brand)	KERN
EAN	4058072136901
Transport temperature	Room temperature
Color	light gray/dark gray
Readability	1 g
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Platform scale
Weighing range MAX	150 kg

