



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN platform scale DS, 36 kg

€1,070.00
plus VAT &
Shipping

Product Images



Description

High-load precision scale, optimal for large-volume or heavy samples that must be weighed with high accuracy
Operation Step by Step via Yes/No dialog in the display Numerical tare deduction for known container weight. Useful for fill quantity checks
Accurate counting: Automatic reference optimization gradually improves the average value of the part weight
Freely programmable weighing unit, e.g., display directly in thread length g/m, paper weight g/m², or similar. Protective work cover included in scope of delivery

Additional Information

No.	KP-4019
Manufacturer (Brand)	KERN
EAN	4058072136963
Transport temperature	Room temperature
Color	light gray dark gray
Readability	200 mg
Operation via	Membrane keypad
Energy supply	Power grid
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Platform scale
Weighing range MAX	36 kg

