



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN platform scale IFS, 150 kg | 300 kg

**Price on
request
plus VAT &
Shipping**

Product Images





Description

- Heavy-duty industrial standard suitable for rough industrial use
- Ergonomic evaluation unit with large keypad and high-contrast LCD displays for convenient input and reading of e.g. tare values, reference weights, limit values etc.
- Accurate counting: Manual reference optimization gradually improves the average value of the part weight
- Totalizing of counted parts
- Three displays for weight indication (verifiable), reference weight, total piece count
- 100 article memory locations for master data such as tare weight (container weight), reference quantity, reference weight, etc.
- Printout with date and time
- Dust cover included in delivery

Additional Information

No.	KP-4467
Manufacturer (Brand)	KERN
EAN	4058072140861
Transport temperature	Room temperature
Color	light gray dark gray
Readability	2 g
Operation via	Membrane keypad
Energy supply	Power grid
Material weighing plate	Stainless steel
Weighing range MAX	300 kg
Manufacturer	CORE
Type scale	Platform scale

