



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale 440, 200 g

**Price on
request
plus VAT &
Shipping**

Product Images



Description

- To make it easier for you to make the right choice for your KERN scale, we have listed a quality code for each model, which is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

Additional Information

No.	KP-0190
Manufacturer (Brand)	KERN
VGKL number	645844033
EAN	4058072124908
Transport temperature	Room temperature
Color	light greylgrey
Wide	165 mm
Height	230 mm
Depth	230 mm
Readability	10 mg
Operation via	Membrane keypad
Weight	1.36 kg
Reproducibility	10 mg
Scales features	with switchable weighing units
Manufacturer	CORE
Material weighing plate	Stainless steel
Type scale	Precision scale
Weighing range MAX	200 g

