




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale 572, 1,600 kg

**€735.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Many typical laboratory functions, such as formulation function, percentage determination, GLP logging and the high precision make the KERN 572 a reliable partner for daily work in the laboratory
- The rugged design, the typical industrial functions like piece counting, vibration-free weighing and the large weighing ranges predestine this series for all industrial applications where highest precision is required
- Freely programmable weighing unit, e.g. display directly in yarn length g/m, paper weight g/m<sup>2</sup>, etc.
- Rugged die-cast aluminum housing provides stability, protects the weighing technology and is insensitive to daily use
- Ring-shaped draft shield standard for models with weighing plate size Ø 106 mm. Weighing chamber Ø×H 157×43 mm
- Eyelet and hook for below-balance weighing standard on models with readability = 0.001 g
- Dust cover included in delivery
- Models with resolution  $\geq 240,000$  points: Level indicator for exact leveling of the balance as standard equipment
- To make it easier for you to make the right choice of your KERN balance, we have listed a quality code for each model. This code consists of two quality features and provides you with a further decision-making aid for the perfect balance for your application in addition to the product features, technical data and pictograms.

## Additional Information

---

|                         |                                |
|-------------------------|--------------------------------|
| No.                     | KP-0653                        |
| Manufacturer (Brand)    | KERN                           |
| VGKL number             | 645857233                      |
| EAN                     | 4058072126315                  |
| Transport temperature   | Room temperature               |
| Color                   | light greylgrey                |
| Wide                    | 180 mm                         |
| Height                  | 310 mm                         |
| Depth                   | 310 mm                         |
| Readability             | 0,01 g                         |
| Operation via           | Membrane keypad                |
| Weight                  | 4.4 kg                         |
| Reproducibility         | 10 mg                          |
| Scales features         | with switchable weighing units |
| Material weighing plate | Stainless steel                |
| Weighing range MAX      | 1 kg                           |
| Manufacturer            | CORE                           |
| Type scale              | Analytical balance             |

