



## KERN precision scale EHA 10001

### Product Images

---



### Description

---

Thanks to the compact and robust design, the bright display as well as the high precision, this series is ideal for use in laboratories, quality control, production as well as school and teaching e.g. biology, chemistry, physics lessons

- Large, shock-resistant stainless steel weighing plate, easy and hygienic to clean
- Particularly flat design
- Ergonomically optimized control panel
- Simple and convenient 5-key operation
- Large, high-contrast display, facilitates reading
- Safe and non-slip installation due to rubber feet
- Level indicator and foot screws for precise leveling of the balance as standard, ensuring the most accurate weighing results
- CAL adjustment program for setting the accuracy, external test weights at extra cost

# Additional Information

No.	KP-9002
Manufacturer (Brand)	KERN
EAN	4058072285685
Color	hellgrau grau
Readability	0.1 g
Operation via	Folientastatur
Reproducibility	0.3 g
Scales features	ultraflach
Material weighing plate	Edelstahl
Calibration type	mit externen Gewichten
Weighing range MAX	1000 g
Manufacturer	KERN
Dimensions weighing plate	120 x 120 mm
Type scale	Präzisionswaage
Weighing units	g

