




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale EHA 10001

€75.00
plus VAT &
Shipping

Product Images



Description

Thanks to the compact and robust design, the bright display as well as the high precision, this series is ideal for use in laboratories, quality control, production as well as school and teaching e.g. biology, chemistry, physics lessons

- Large, shock-resistant stainless steel weighing plate, easy and hygienic to clean
- Particularly flat design
- Ergonomically optimized control panel
- Simple and convenient 5-key operation
- Large, high-contrast display, facilitates reading
- Safe and non-slip installation due to rubber feet
- Level indicator and foot screws for precise leveling of the balance as standard, ensuring the most accurate weighing results
- CAL adjustment program for setting the accuracy, external test weights at extra cost

Additional Information

No.	KP-9002
Manufacturer (Brand)	KERN
EAN	4058072285685
Transport temperature	Room temperature
Color	light grey grey
Scales features	ultra flat
Readability	0,1 g
Operation via	Membrane keypad
Manufacturer	CORE
Material weighing plate	Stainless steel
Dimensions weighing plate	120 x 120 mm
Reproducibility	0,3 g
Calibration type	with external weights
Type scale	Precision scale
Weighing range MAX	1000 g
Weighing units	g

