



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale EHA 30000

**€70.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Thanks to the compact and robust design, the bright display as well as the high precision, this series is ideal for use in laboratories, quality control, production as well as school and teaching e.g. biology, chemistry, physics lessons

- Large, shock-resistant stainless steel weighing plate, easy and hygienic to clean
- Particularly flat design
- Ergonomically optimized control panel
- Simple and convenient 5-key operation
- Large, high-contrast display, facilitates reading
- Safe and non-slip installation due to rubber feet
- Level indicator and foot screws for exact leveling of the balance as standard, ensuring the most accurate weighing results
- CAL adjustment program for setting the accuracy, external test weights at extra cost

## Additional Information

---

No.	KP-9003
Manufacturer (Brand)	KERN
EAN	4058072285692
Transport temperature	Room temperature
Color	light greylgrey
Readability	0,1 g
Operation via	Membrane keypad
Reproducibility	0,3 g
Scales features	ultra flat
Material weighing plate	Stainless steel
Calibration type	with external weights
Weighing range MAX	3000 g
Manufacturer	CORE
Dimensions weighing plate	120 x 120 mm
Type scale	Precision scale
Weighing units	g

