



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale EW , 620 g

€1,080.00
plus VAT &
Shipping

Product Images



Description

- KERN EG-N with internal automatic adjustment in case of temperature changes or time-controlled according to defined intervals. Guarantees high accuracy and makes independent of location
- eW-NM only: CAL adjustment program for setting the accuracy by means of an external test weight
- Stable temperature behavior
- Short settling time
- High mechanical robustness
- High corner load safety
- Capacity display: A rising luminous band indicates the still available weighing capacity
- GLP/ISO logging of weighing values
- Totalizing of counting parts
- Draft shield included as standard for models with weighing plate size Ø 118 mm, weighing chamber W×D×H 158×130×78 mm
- Dust cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code is made up of two quality features and, in addition to the product features, technical data and pictograms, provides you with a further decision-making aid for the perfect scale for your application.

Additional Information

| | |
|-------------------------|--------------------------------|
| No. | KP-7315 |
| Manufacturer (Brand) | KERN |
| VGKL number | 645816203 |
| EAN | 4058072284213 |
| Transport temperature | Room temperature |
| Color | light greylgrey |
| Operation via | Membrane keypad |
| Weight | 2.96 kg |
| Reproducibility | 1 mg |
| Scales features | with switchable weighing units |
| Material weighing plate | Stainless steel |
| Weighing range MAX | 620 g |
| Manufacturer | CORE |
| Type scale | Precision scale |

