




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale EWJ, 300 g

**€600.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Convenient operating philosophy: all main functions have their own key on the control panel
- Internal automatic calibration timed every 2 hours. Guarantees high accuracy and makes the scale independent of location
- Capacity indicator: A rising luminous band shows the remaining available weighing capacity
- Draft shield standard on models with = 600 g weighing capacity W×D×H 134×128×80 mm
- eWJ 300-3H only: Large glass draft shield with 3 sliding doors for convenient access to the sample, weighing chamber W×D×H 155×175×217 mm
- RS-232 data interface as standard (not for EWJ-SM)
- USB data interface for transferring weighing data to PC, printer, USB sticks etc. (not for EWJ-SM)
- Dust cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect scale for your application.

## Additional Information

---

|                         |                                |
|-------------------------|--------------------------------|
| No.                     | KP-4155                        |
| Manufacturer (Brand)    | KERN                           |
| VGKL number             | 105839303                      |
| EAN                     | 4058072138059                  |
| Transport temperature   | Room temperature               |
| Color                   | light greylgrey                |
| Readability             | 1 mg                           |
| Operation via           | Membrane keypad                |
| Weight                  | 3.8 kg                         |
| Reproducibility         | 2 mg                           |
| Scales features         | with switchable weighing units |
| Manufacturer            | CORE                           |
| Material weighing plate | Stainless steel                |
| Type scale              | Precision scale                |
| Weighing range MAX      | 300 g                          |

