



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale EWJ, 6 kg

€480.00
plus VAT &
Shipping

Product Images



Description

- Convenient operating philosophy: all main functions have their own key on the control panel
- Internal automatic calibration timed every 2 hours. Guarantees high accuracy and makes the scale independent of location
- Capacity indicator: A rising luminous band shows the remaining available weighing capacity
- Draft shield standard on models with = 600 g weighing capacity W×D×H 134×128×80 mm
- eWJ 300-3H only: Large glass draft shield with 3 sliding doors for convenient access to the sample, weighing chamber W×D×H 155×175×217 mm
- RS-232 data interface as standard (not for EWJ-SM)
- USB data interface for transferring weighing data to PC, printer, USB sticks etc. (not for EWJ-SM)
- Dust cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect scale for your application.

Additional Information

No.	KP-4158
Manufacturer (Brand)	KERN
VGKL number	105839305
EAN	4058072138080
Transport temperature	Room temperature
Color	light greylgrey
Readability	100 mg
Operation via	Membrane keypad
Weight	4.48 kg
Reproducibility	100 mg
Material weighing plate	Stainless steel
Weighing range MAX	6 kg
Manufacturer	CORE
Type scale	Precision scale
Weighing units	g

