




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN precision scale EWJ, 6 kg

**€480.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- Convenient operating philosophy: all main functions have their own key on the control panel
- Internal automatic calibration timed every 2 hours. Guarantees high accuracy and makes the scale independent of location
- Capacity indicator: A rising luminous band shows the remaining available weighing capacity
- Draft shield standard on models with = 600 g weighing capacity W×D×H 134×128×80 mm
- eWJ 300-3H only: Large glass draft shield with 3 sliding doors for convenient access to the sample, weighing chamber W×D×H 155×175×217 mm
- RS-232 data interface as standard (not for EWJ-SM)
- USB data interface for transferring weighing data to PC, printer, USB sticks etc. (not for EWJ-SM)
- Dust cover included in delivery
- To make it easier for you to make the right choice of your KERN scale, we have listed a quality code for each model. This code consists of two quality features and provides you, in addition to the product features, technical data and pictograms, with a further decision-making aid for the perfect scale for your application.

## Additional Information

---

No.	KP-4158
Manufacturer (Brand)	KERN
VGKL number	105839305
EAN	4058072138080
Transport temperature	Room temperature
Color	light greylgrey
Readability	100 mg
Operation via	Membrane keypad
Weight	4.48 kg
Manufacturer	CORE
Material weighing plate	Stainless steel
Reproducibility	100 mg
Type scale	Precision scale
Weighing range MAX	6 kg
Weighing units	g

