



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN precision scale KB, 10 kg

Price on
request
plus VAT &
Shipping

Product Images



Description

Many typical laboratory functions, such as formulation function, percentage determination, GLP logging and the high precision make the KERN KB a reliable partner for daily work in the laboratory Operation step by step via yes/no dialog in the display Freely programmable weighing unit, e.g. display directly in thread length g/m, paper weight g/m², or similar. Ring-shaped draft shield as standard, only for models with weighing plate size Ø 81 mm, weighing chamber Ø×H 90×40 mm Dust cover included Models with resolution >= 240,000 points and models with verification approval: Level indicator for exact leveling of the balance as standard

Additional Information

No.	KP-4539
Manufacturer (Brand)	KERN
VGKL number	105838105
EAN	4058072141493
Transport temperature	Room temperature
Color	light greylgrey
Readability	50 mg
Operation via	Membrane keypad
Weight	2.2 kg
Reproducibility	50 mg
Scales features	with switchable weighing units
Material weighing plate	Stainless steel
Weighing range MAX	10 kg
Manufacturer	CORE
Type scale	Precision scale

